

Proximel[®]

Care. Closest to skin.




Proximel[®] helps you manage
patient wounds easily, efficiently
and with confidence.


Whether used for prevention of pressure injuries or wound care, Proximel® is gentle on patients.

Proximel® silicone wound dressings feature a unique five-layer design that promotes healing and supports pressure injury prevention by quickly and effectively managing wound exudate while reducing and redistributing pressure at the injury site. Part of the HARTMANN family of wound care solutions, Proximel® is a versatile choice for a wide variety of care settings, prevention protocols and wound types.


A five-layer design that improves upon traditional silicone dressings




Proximel's distinctive properties allow you to confidently **lift** and **reposition** the dressing while reducing the risk of adherence failure.




High MVTR rate supports high-absorption capacity and good moisture balance.



Superabsorbent core supports **fast absorption** and helps minimize maceration risk.



Reduces and redistributes pressure across the surface to reduce both peak and overall level of force.



A versatile solution

Proxigel® is ideal for use throughout various care settings for both the prevention and treatment of wounds.⁵



Proxigel®:
A soft, five-layer
silicone bordered
dressing



Ordering Information

Proxigel®

Order Number	Description	Dressings per Box
14100000	Proxigel Silicone Dressing with Border, 3" x 3"	10
14200000	Proxigel Silicone Dressing with Border, 4" x 4"	10
14300000	Proxigel Silicone Dressing with Border, 5" x 5"	10
14400000	Proxigel Silicone Dressing with Border, 6" x 6"	5
14500000	Proxigel Silicone Dressing with Border, 4" x 12"	5
14600000	Proxigel Silicone Dressing with Border, Small Sacrum	5
14700000	Proxigel Silicone Dressing with Border, Large Sacrum	5
14800000	Proxigel Silicone Dressing with Border, Heel	5

Contact your local HARTMANN representative,
or visit us at www.hartmannusa.com for more information.

HARTMANN USA, Inc.
Rock Hill, SC 29730
1-800-243-2294

XLIT 2757 Rev. 2 (0118)

1. Internal lab data test conducted with Proxigel, Allevyn Gentle Border and Mepilex Border.
2. Testing performed per BS EN 13726-2:2002.
3. Testing performed per BS EN 13726-2:2002.
4. Simulated with steel ball, data on file.
5. Emily Haesler, editor. National Pressure Ulcer Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers: Quick Reference Guide. Perth: Cambridge Media, 2014.



Going further
for health